

Mite Extract Dpb

1. Product

Appearance :	beige-colored powder
Origin :	<i>Dermatophagoides pteronyssinus</i>
Preparation :	Extract from mite bodies, freeze-dried
Lot :	
Volume :	2 mg (total protein, BCA assay)

2. Storage

store at -20 °C

3. Notes

For research use only. Not for any clinical, therapeutic, or diagnostic use.

4. Manufacture & Retail

Biostir Inc.



Annex Shinosaka 4A, 2-9-1 Miyahara, Yodogawa-ku, Osaka, Osaka, Japan
tel +81-6-6151-3950 fax +81-6-6151-3940

e-mail: info@biostir.com

URL: <https://www.biostir.com/>